# PROXIMATE® Skin Staplers

# The number one choice in skin staplers\*

The **PROXIMATE® Rotating Head** Skin Stapler features a 360° rotating head assembly for improved visibility and a ratchet mechanism for easier staple placement.

The **PROXIMATE® PX Fixed Head** Skin Stapler provides many of the same features as the PROXIMATE Rotating Head Stapler, but in a fixed-head format.

The **PROXIMATE® Plus MD** Skin Stapler is a high-value, low-cost skin stapler that permits multi-directional release in an ergonomic design.



PROXIMATE Rotating Head Stapler

## PROXIMATE Rotating Head Stapler

- Head rotates 360° for improved visibility and access
- Ratchet mechanism for easy staple placement
- Strong, rectangular stainless steel staples resist opening and minimize staple rotation
- Ergonomic pistol-grip handle is intuitive and comfortable to use

## PROXIMATE PX Fixed Head Stapler

- Ratchet Mechanism for easy staple placement
- Ergonomic pistol-grip handle is intuitive and comfortable to use

# PROXIMATE Plus MD Stapler

- Improved kick-off spring design for multi-directional release
- Ergonomic pistol-grip handle is intuitive and comfortable to use
- Alignment indicator improves visibility



PROXIMATE PX Fixed Head Stapler



PROXIMATE Plus MD Stapler

# Ordering Information

## **Skin Staplers**

Description	Device	Device Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Qty/box
PROXIMATE® Rotating Head Skin Stapler	O 2000001 to 36	PRR35	35	3.9mm	0.53mm	5.7mm	6
		PRW35	35	3.9mm	0.58mm	6.9mm	6
PROXIMATE® PX Fixed Head Skin Stapler		PXR35	35	3.9mm	0.53mm	5.7mm	6
		PXW35	35	3.9mm	0.58mm	6.9mm	6
PROXIMATE® Plus MD Skin Stapler	350	PMR35	35	3.9mm	0.53mm	5.7mm	6
		PMW35	35	3.9mm	0.58mm	6.9mm	6
Description	Device	Device Code	Qty/box				
Squeeze-Handle Staple Extractor		PSX	12				

Squeeze-Handle Staple Extractor is compatible with all PROXIMATE Skin Staplers

## How to Order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

#### Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

## Johnson & Johnson Gateway: www.jnjgateway.com/commerce

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

### Global Healthcare Exchange (GHX): www.ghx.com

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

Electronic Data Interchange (EDI): JJHCS EDI Help Line: 1-800-262-2888

#### Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

# **Customer Support**

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at **customersupport@eesus.jnj.com** or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.