

# Intersurgical EcoLite™ adult medium concentration oxygen mask

The lightest and most comfortable mask available







Oxygen and Aerosol Therapy • Oxygen Therapy







## Intersurgical ECOLite<sup>™</sup> adult medium concentration oxygen mask

#### Welcome to the Comfort Zone

Following the outstanding success of the original Intersurgical EcoLite adult oxygen mask, with millions sold around the world, we are pleased to announce the arrival of the next generation.

Incorporating new features such as an extended chin with two integral seals, the new mask ensures a comfortable fit to a wider range of face shapes and sizes.

The world's lightest and most comfortable oxygen mask has just become even more comfortable.

#### Comfortable for the patient

Incurved nose seal

No metal nose clip

**Exhalation ports** 

Designed to minimise the amount of rebreathed carbon dioxide

MRI compatible

Conforms to different nose

shapes designed to prevent oxygen entering patient's eyes

The unique features of the original Intersurgical EcoLite oxygen mask have been retained and include the following:

#### Comfortable for the environment

Various groups have questioned the use of PVC in medical products and it's impact on the environment. Intersurgical's aim is to reduce the environmental impact of its products and processes. The utilisation of new materials and the latest manufacturing technology has resulted in the elimination of PVC from the mask shell resulting in a reduced environmental impact.





www.intersurgical.com/info/eco

Video available at www.intersurgical.com

#### Lighter

50% lighter than a PVC mask, more comfortable for the patient

#### Greener @

The elimination of PVC from the masks has dramatically lowered their environmental impact

#### Soft face seals

Increased patient comfort

### Elastic can be positioned under or over the ears


Below ear position eliminates trauma to top of ears

"under chin"

position

Large

chin seal





#### A choice of "under chin" or "on chin" positions

Provides a better fit on a wider range of patient face shapes



Code Description Box Qtv. 1136015 Intersurgical EcoLite™, adult, medium concentration oxygen mask

"on chin"

position

Small

chin seal

1135015 Intersurgical EcoLite™, adult, medium concentration oxygen mask with tube, 2.1m





Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK T: +44 (0)118 9656 300 F: +44 (0)118 9656 356 info@intersurgical.com www.intersurgical.com

UK • Ireland • France • Germany • Spain • Portugal Italy • Benelux • Sweden • Denmark • Lithuania • Russia Czech Republic • Turkev • South Africa • China • Japan Taiwan • Philippines • USA • Canada • Colombia • Australia







